



## SAMPLING PROCEDURE MANUAL

WRITTEN BY: DR. METİN İLEMİN

## SAMPLING PROCEDURE MANUAL

### 1. INTRODUCTION AND CLINICAL HISTORY


Diagnosing a fish disease requires more than just observing clinical signs. **Many diseases share common symptoms** (e.g., darkening, exophthalmia/pop-eyes, ascites, hemorrhage, splenomegaly, changes in organ color/shape). A complete diagnosis relies on combination of physical examinations, sampling, and environmental data.

**Key Clues for Diagnosis:** Observation of fish behavior, reflexes, lesion description, necropsy, bacteriology, and parasitological examination.

**Fundamental Context:** Record information on:

- **Transportation** and **handling** history.
- **Previous treatments** (medicated fish should be avoided if possible).
- Possible **stressful events** (heavy rain, predators, power outages, etc.).
- Current **water parameters**.

### 2. SAMPLE SELECTION: CHOOSING THE RIGHT FISH




The quality of the sample is must for an accurate diagnosis.

**Select:** **Sick fish** showing any disease symptom or **moribund fish** (those near death).

- Sick fish are often easy to spot: they may be lethargic, separated from the group, or found near the sides, surface, or inlet/outlet pipes.
- **Avoid: Never sample dead fish.** Tissue condition, chemistry, and ectoparasites rapidly deteriorate or detach post-mortem, leading to a false picture of the disease state.

### 3. BACTERIOLOGY SAMPLING

**NOTE: Aseptic procedures** must be used for all bacteriology samples to prevent contamination. Samples for bacteriology must be taken **before** any other samples (e.g., PCR, histopathology) to avoid cross-contamination from one organ to another.



#### Procedure

- 1. Preparation:** Place the selected fish on its **right flank**.
- 2. Incision:** Using a **sterile scalpel** (new or cleaned with alcohol/bleach), make an incision in the **abdominal area**. **Carefully avoid perforating the digestive system**.
- 3. Extending the Cut:** Use a pair of **sterile scissors** to extend the incision:
  - Forward to the **operculum**.
  - Backward to the **anus**.
- 4. Exposing Organs:** Using the sterile scissors, **carefully cut away the musculature** overlaying the internal organs to expose the intestine and the underlying organs (kidney and spleen).
- 5. Swabbing the Organs:** Using a **sterile bacterial swab**, gently push in and insert the tip into the **kidney**.
- 6. Repeat:** Perform the **same swabbing procedure** for the **spleen**.
- 7. Packaging:** **Close the tap/cap of the sterile swab** securely and prepare the sample for **immediate shipment** to the diagnostic laboratory following their instructions.

#### 4. FTA CARD SAMPLING

FTA cards capture and preserve DNA/RNA directly in the field by binding cells onto a chemically treated fiber matrix. The card instantly lyses and inactivates pathogens while stabilizing nucleic acids at room temperature, eliminating the need for cold-chain transport. This makes sampling easy, safe, and reliable even for non-expert farm staff. Once in the lab, a small piece from the card is enough for high-quality PCR analysis, providing a fast, practical, and contamination-resistant solution for aquaculture diagnostics.



#### 5. FTA CARD USAGE AND CAPACITY

FTA cards are specialized paper that **absorbs and fixes genetic material** (nucleic acids) to the cellulose upon contact with tissue. The nucleic acids are immobilized and preserved, allowing for sample transport at **room temperature**.

- Each FTA card contains **4 circles or wells**.
- Up to **3 samples** can be pooled in each well.
- Total capacity: **12 samples** per card (4 wells 3 pooled samples).  
Pooling allows multiple fish to be treated as a single sample:
- **Small Fish (Fry or Fingerlings)**: 5 to 10 fish may be pooled together
- **Larger Fish**: Tissues from a maximum of **5 fish** may be pooled.

Prior to shipment, **always complete the sample submission form** and notify the recipient laboratory of the impending submission, following their advice regarding labeling and shipping dates.

## 6. SAMPLING METHOD

Use the same procedure described above for bacteriology. Collect the sample using a foam applicator or a sterile swab, then apply the material from the internal organs onto the FTA card. Alternatively, you may directly excise a small piece of the organ and gently imprint it onto the FTA Card.



## 7. SHIPMENT

FTA cards make sample transport and shipment significantly easier and simpler due to several key features:

- **No Biohazard Classification:** The cards actively **inactivate the pathogen**, ensuring the sample is **not categorized as a Biological Sample Type B**. This crucial feature eliminates the need for extensive paperwork, such as import licenses.
- **Temperature Stability:** Samples can be transported at **ambient (room) temperature**.
- **Extended Shelf Life:** The shipment window is long, allowing transport for **up to 4 weeks**.
- **Simple Mailing:** You can use **normal mail** for shipment by simply packing the cards securely in an envelope.